

Bead-Beater DNA Preparation Protocol

1. Pellet cells/material for 5 minutes at 10,000g
2. Resuspend pellet in 800 μ l lysis solution
3. Add 0.8g of 0.1mm glass beads (or equal volume)
4. Disrupt for 30 seconds at maximum speed.
 - a. Time many vary depending upon starting material/organisms present
5. Centrifuge at 10,000g for 3 minutes
6. Transfer supernatant to new 1.5ml eppendorf tube.
7. Add equal volume phenol/chloroform.
8. Incubate at 68°C for 5 minutes with 2-3 mixes.
9. Centrifuge at 10,000g for 3 minutes.
10. Transfer aqueous phase to new 1.5ml eppendorf.
11. (OPTIONAL) Do another phenol/chloroform extraction.
 - a. Dependent upon previous 260/280 readings from similar sample types.
12. Add 2 volumes 95% EtOH and 1/10th volume NaAcetate (3M)
13. Precipitate at -20°C O/N.
14. Centrifuge at 10,000g for 45-60 minutes at 4°C.
15. Wash with 70% EtOH.
16. Centrifuge 10,000g for 30 minutes at 4°C.
17. Air dry pellet
18. Resuspend in DNase/RNase-free TE (pH 7.0 or 8.0)

Lysis Solution

2.7g NaAcetate (33mM)
5.0g SDS (17.5mM)
0.34g EDTA (~3mM)
pH 5.5
q.s. to 1 Liter total volume